# EU Power Adapter Infrared Therapy Light Wavelength 660nm 850nm 1400W **Upgrade Your Health Routine**



for more
products
please
visit us on g
on grov
vled-ligh
ht.com

Basic Information	

- Brand Name: ANROFU · Certification:
  - FDA/FCC/CE/ROHS/UL/SAA/C TICK
- Model Number: M1400 1
- Minimum Order Quantity:
- L904\*W385\*H130mm 10KG Packaging Details:
- Delivery Time: 5-8 days
- T/T,Western Union • Payment Terms:
- Supply Ability: 1000pcs/ day



# **Product Specification**

- Input Voltage: AC100~240V • Power Adapter: US, EU, AU, UK
- L764\*W310\*H60mm • Dimensions:
- Wavelength: 280LED 660nm\*140pcs+850nm\*140pcs
- Treatment Time: 5-30 Minutes
  - Infrared Light Therapy Device
  - 1400W
- "-20~40 45%~95%RH" • Temperature: 1400W Infrared Therapy Light,
  - 660nm Infrared Therapy Light, 850nm Infrared Therapy Light



# More Images

• Type: • Power:

• Highlight:



# Irradiance & Coverage







# Our Product Introduction

ನ್

Features:

Product Description

Product Name: Infrared Therapy Light Input Voltage: AC100-240V Wavelength: 660nm 850nm Power Output: 1400W Type: Infrared Light Therapy Device Power Supply: AC Adapter

### **Applications:**

The ANROFU M1400 Infrared Therapy Light, originating from China, is a versatile product suitable for various applications and scenarios, thanks to its advanced features and benefits

One of the key uses of this therapy light is for muscle recovery and pain relief. The IR Infrared LED technology with a wavelength of 660nm 850nm penetrates deep into the muscles, providing soothing relief from aches, sprains, and stiffness. Athletes, fitness enthusiasts, or individuals with muscle-related issues can benefit greatly from this Infrared Heat Lamp.

or more products please visit us on growled-light.com

Moreover, the ANROFU M1400 is ideal for skincare treatments. The Infrared Red Light Therapy helps stimulate collagen production, improve skin elasticity, and reduce signs of aging. It can be used in beauty salons, spas, or even at home for rejuvenating skincare routines.

Individuals looking for relaxation and stress relief can also find this therapy light beneficial. The gentle heat emitted by the device can promote relaxation, relieve tension, and create a calming atmosphere. The auto shut-off feature ensures safety and peace of mind during use, making it suitable for meditation or relaxation sessions lasting between 5 to 30 minutes.

With compact dimensions of L764\*W310\*H60MM, the ANROFU M1400 is portable and easy to use in various settings. Whether at home, in a clinic, or while traveling, this therapy light can be powered using an AC/DC adapter, offering convenience and flexibility.

In conclusion, the ANROFU M1400 Infrared Therapy Light is a versatile and effective solution for muscle recovery, skincare treatments, relaxation, and pain relief. Its advanced features, including the Infrared Heat Lamp technology and auto shut-off function, make it a valuable addition to wellness routines in a range of settings.

### **Customization:**

Enhance your infrared light therapy experience with our customizable ANROFU M1400 Infrared Light Therapy Device, originating from CHINA. Designed for optimal performance, this device offers treatment times ranging from 5 to 30 minutes, making it suitable for various needs.

Powered by an AC adapter, this device provides a power output of 1400W, ensuring consistent and reliable operation. The large treatment area of up to 20 inches caters to a wide range of users, allowing you to target specific areas effectively.

Experience the benefits of Infrared Light Therapy Devices with our versatile ANROFU M1400 model, suitable for Infrared LED Light Therapy and Infrared Laser Therapy Machine applications.

### Support and Services:

Our Infrared Therapy Light product comes with comprehensive technical support and services to ensure you have a seamless experience. Our team of experts is available to assist you with any product-related queries or issues you may encounter.

Our services include troubleshooting assistance, product maintenance tips, and guidance on optimal usage of the Infrared Therapy Light. We are committed to providing you with the support you need to make the most out of your product.

For any technical support or service inquiries, please do not hesitate to reach out to our dedicated team who will be more than happy to assist you.

ANROFU Shenzhen Hanyuan Intelligent Technology Co., Ltd.

C +86 15017943686 S wanghaiyan@anrofu.com S growled-light.com

3 Floor No.70, Zhulong Load, Fourth Industrial Shiyan Street Bao'An District Shenzhen City, Guangdong, China